

## VE-Cadherin (3G8) Mouse mAb

Cat No: HR1AM2495

For research use only

### Overview

Product Name	VE-Cadherin (3G8) Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	IHC 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	80-115kDa
GeneID?Human?	1003
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD144,Cadherin 5, VEC,CDH5
Background	In endothelial cells, VE-cadherin signaling, expression, and localization correlate with vascular permeability and tumor angiogenesis.