




CGI-43 Polyclonal Antibody

Cat No: HR1AP6345

For research use only

Overview

Product Name	CGI-43 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GenelID?Human?	CCZ1/C7orf28B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCZ1; C7orf28A; CGI-43; Vacuolar fusion protein CCZ1 homolog; CCZ1B; C7orf28B; Vacuolar fusion protein CCZ1 homolog B; Vacuolar fusion protein CCZ1 homolog-like
Background	similarity:Belongs to the UPF0550 family.,



bioelsa