




CCDC4 Polyclonal Antibody

Cat No: HR1AP2997

For research use only

Overview

Product Name	CCDC4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GenID?Human?	BEND4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BEND4; CCDC4; BEN domain-containing protein 4; Coiled-coil domain-containing protein 4
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 BEN domain.,



bioelsa