

CCDC45 Polyclonal Antibody

Cat No: HR1AP2998

For research use only

Overview

Product Name	CCDC45 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	95kDa
GeneID?Human?	CEP95
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEP95; CCDC45; CEP45; Centrosomal protein of 95 kDa; Coiled-coil domain-containing protein 45
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,