



BUD31 Polyclonal Antibody

Cat No: HR1AP2913

For research use only

Overview

Product Name	BUD31 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	17kDa
GenID?Human?	BUD31
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BUD31; EDG2; Protein BUD31 homolog; Protein EDG-2; Protein G10 homolog
Background	similarity:Belongs to the BUD31 (G10) family.,

