



Brp16 Polyclonal Antibody

Cat No: HR1AP2903

For research use only

Overview

Product Name	Brp16 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GenelD?Human?	FAM203A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FAM203A; BRP16; C8orf30A; Protein FAM203A; Brain protein 16
Background	similarity:Belongs to the UPF0570 (BRP16/HGH1) family.,