

AARSD1 Polyclonal Antibody

Cat No: HR1AP2669

For research use only

Overview

Product Name	AARSD1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	AARSD1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	AARSD1; Alanyl-tRNA editing protein Aarsd1; Alanyl-tRNA synthetase domain-containing protein 1
Background	caution:The sequence shown here is derived from an Ensembl automatic analysis pipeline and should be considered as preliminary data.,similarity:Belongs to the class-II aminoacyl-tRNA synthetase family. AARSD1 subfamily.,