



## ALY Polyclonal Antibody

Cat No: HR1AP5547

For research use only

### Overview

Product Name	ALY Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	27kDa
GenID?Human?	ALYREF
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ALYREF; ALY; BEF; THOC4; THO complex subunit 4; Tho4; Aly of AML-1 and LEF-1; Aly/REF export factor; Transcriptional coactivator Aly/REF; bZIP-enhancing factor BEF
Background	Aly/REF export factor(ALYREF) Homo sapiens The protein encoded by this gene is a heat stable, nuclear protein and functions as a molecular chaperone. It is thought to regulate dimerization, DNA binding, and transcriptional activity of basic region-leucine zipper (bZIP) proteins. [provided by RefSeq, Jul 2008],