

C/EBP ? (phospho Ser21) Polyclonal Antibody

Cat No: HR1AP2579

For research use only

Overview

Product Name	C/EBP ? (phospho Ser21) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	CEBPA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CEBPA; CCAAT/enhancer-binding protein alpha; C/EBP alpha
Background	CCAAT/enhancer binding protein alpha(CEBPA) Homo sapiens This intronless gene encodes a transcription factor that contains a basic leucine zipper (bZIP) domain and recognizes the CCAAT motif in the promoters of target genes. The encoded protein functions in homodimers and also heterodimers with CCAAT/enhancer-binding proteins beta and gamma. Activity of this protein can modulate the expression of genes involved in cell cycle regulation as well as in body weight homeostasis. Mutation of this gene is associated with acute myeloid leukemia. The use of alternative in-frame non-AUG (GUG) and AUG start codons results in protein isoforms with different lengths. Differential translation initiation is mediated by an out-of-frame, upstream open reading frame which is located between the GUG and the first AUG start codons. [provided by RefSeq, Dec 2013],