



## LDHD Polyclonal Antibody

Cat No: HR1AP9264

For research use only

### Overview

Product Name	LDHD Polyclonal Antibody
Source	Rabbit
Applications	WB
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	55kDa
GenID?Human?	LDHD
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Probable D-lactate dehydrogenase, mitochondrial (DLD) (Lactate dehydrogenase D) (EC 1.1.2.4)
Background	lactate dehydrogenase D(LDHD) Homo sapiens The protein encoded by this gene belongs to the D-isomer specific 2-hydroxyacid dehydrogenase family. The similar protein in yeast has both D-lactate and D-glycerate dehydrogenase activities. Alternative splicing occurs at this locus and two transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Jul 2008],