

RRP7A Polyclonal Antibody

Cat No: HR1AP4954

For research use only

Overview

Product Name	RRP7A Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	RRP7A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RRP7A; CGI-96; Ribosomal RNA-processing protein 7 homolog A; Gastric cancer antigen Zg14; RRP7B; Putative ribosomal RNA-processing protein 7 homolog B; Putative gastric cancer antigen Zg14-like protei
Background	caution:Could be the product of a pseudogene.,similarity:Belongs to the RRP7 family.,similarity:Contains 1 RRM (RNA recognition motif) domain.,