

# TFIIIC102 Polyclonal Antibody

**Cat No: HR1AP4700**

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | TFIIIC102 Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                         |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                |
| Applications          | WB,IHC-p,ELISA                                                                                                                                                                                                                                                                                                                                                                                                        |
| Species Reactivity    | Human,Mouse                                                                                                                                                                                                                                                                                                                                                                                                           |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                          |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                               |
| Observed band         | 101kDa                                                                                                                                                                                                                                                                                                                                                                                                                |
| GeneID?Human?         | GTF3C3                                                                                                                                                                                                                                                                                                                                                                                                                |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                       |
| Alternative Names     | GTF3C3; General transcription factor 3C polypeptide 3; Transcription factor IIIC 102 kDa subunit; TFIIIC 102 kDa subunit; TFIIIC102; Transcription factor IIIC subunit gamma; TF3C-gamma                                                                                                                                                                                                                              |
| Background            | general transcription factor IIIC subunit 3(GTF3C3) Homo sapiens The protein encoded by this gene is part of the TFIIIC2 complex, which binds to the promoters of small nuclear and cytoplasmic RNA genes in order to recruit RNA polymerase III. The TFIIIC2 complex is composed of six subunits. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011], |