



TFIIC102 Polyclonal Antibody

Cat No: HR1AP4700

For research use only

Overview

Product Name	TFIIC102 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	101kDa
GenelID?Human?	GTF3C3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GTF3C3; General transcription factor 3C polypeptide 3; Transcription factor IIIC 102 kDa subunit; TFIIC 102 kDa subunit; TFIIC102; Transcription factor IIIC subunit gamma; TF3C-gamma
Background	general transcription factor IIIC subunit 3(GTF3C3) Homo sapiens The protein encoded by this gene is part of the TFIIC2 complex, which binds to the promoters of small nuclear and cytoplasmic RNA genes in order to recruit RNA polymerase III. The TFIIC2 complex is composed of six subunits. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Oct 2011],