

TAP-Tag Mouse mAb

Cat No: HR1AM2238

For research use only

Overview

Product Name	TAP-Tag Mouse mAb
Source	Mouse
Applications	WB, IP
Species Reactivity	N/A
Recommended Dilutions	WB 1:5,000-10,000 IP 1:200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	N/A
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	tandem affinity purification
Background	Epitope mapping within the calmodulin binding peptide (CBP) region of the TAP (tandem affinity purification) tag of human.