

RNF113B Polyclonal Antibody

Cat No: HR1AP4496

For research use only

Overview

Product Name	RNF113B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	RNF113B
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RNF113B; RNF161; ZNF183L1; RING finger protein 113B; Zinc finger protein 183-like 1
Background	similarity:Contains 1 C3H1-type zinc finger.,similarity:Contains 1 RING-type zinc finger.,