

Collagen III Mouse mAb Cat No: HR1AM2234

For research use only

Overview

Product Name	Collagen III Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:1,000 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	138kDakDa
GeneID?Human?	1281
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CO3A1 antibody, Collagen alpha 1(III) chain antibody, Collagen III alpha 1 chain precursor antibody, Collagen type III alpha 1 antibody
Background	Collagen alpha-1(III) chain is a fibrillar collagen that is found in extensible connective tissues such as skin, lung, and the vascular system, frequently in association with type I collagen.