



CALHM1 Polyclonal Antibody

Cat No: HR1AP9329

For research use only

Overview

Product Name	CALHM1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	
Observed band	40kDa
GeneID?Human?	CALHM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Calcium homeostasis modulator protein 1; Protein FAM26C
Background	calcium homeostasis modulator 1(CALHM1) Homo sapiens This gene encodes a calcium channel that plays a role in processing of amyloid-beta precursor protein. A polymorphism at this locus has been reported to be associated with susceptibility to late-onset Alzheimer's disease in some populations, but the pathogenicity of this polymorphism is unclear.[provided by RefSeq, Mar 2010],