

Kv10.2 Polyclonal Antibody

Cat No: HR1AP9452

For research use only

Overview

Product Name	Kv10.2 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	KCNH5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Potassium voltage-gated channel subfamily H member 5 (Ether-a-go-go potassium channel 2) (hEAG2) (Voltage-gated potassium channel subunit Kv10.2)
Background	potassium voltage-gated channel subfamily H member 5(KCNH5) Homo sapiens This gene encodes a member of voltage-gated potassium channels. Members of this family have diverse functions, including regulating neurotransmitter and hormone release, cardiac function, and cell volume. This protein is an outward-rectifying, noninactivating channel. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jul 2013],