



Mash2 Polyclonal Antibody

Cat No: HR1AP7367

For research use only

Overview

Product Name	Mash2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	25kDa
GeneID?Human?	ASCL2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ASCL2; BHLHA45; HASH2; Achaete-scute homolog 2; ASH-2; hASH2; Class A basic helix-loop-helix protein 45; bHLHa45; Mash2
Background	achaete-scute family bHLH transcription factor 2(ASCL2) Homo sapiens This gene is a member of the basic helix-loop-helix (BHLH) family of transcription factors. It activates transcription by binding to the E box (5'-CANNTG-3'). Dimerization with other BHLH proteins is required for efficient DNA binding. Involved in the determination of the neuronal precursors in the peripheral nervous system and the central nervous system. [provided by RefSeq, Jul 2008],