

BTBD6 Polyclonal Antibody

Cat No: HR1AP9042

For research use only

Overview

Product Name	BTBD6 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	BTBD6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	BTBD6; BDPL; BTB/POZ domain-containing protein 6; Lens BTB domain protein
Background	similarity:Contains 1 BTB (POZ) domain.,tissue specificity:Expressed in lens.,