



## ECRG4 Polyclonal Antibody

Cat No: HR1AP8868

For research use only

### Overview

Product Name	ECRG4 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p,IF,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GenID?Human?	C2orf40
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	C2orf40; ECRG4; Augurin; Esophageal cancer-related gene 4 protein
Background	function:Probable hormone.,similarity:Belongs to the augurin family.,