



Rrn3 (phospho Ser649) Polyclonal Antibody

Cat No: HR1AP7782

For research use only

Overview

Product Name	Rrn3 (phospho Ser649) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	70kDa
GeneID?Human?	RRN3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RRN3; TIFIA; RNA polymerase I-specific transcription initiation factor RRN3; Transcription initiation factor IA; TIF-IA
Background	function:Required for efficient transcription initiation by RNA polymerase I.,similarity:Belongs to the RRN3 family.,