

CD25 Mouse mAb

Cat No: HR1AM2218

For research use only

Overview

Product Name	CD25 Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human
Recommended Dilutions	IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	3559
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IL2 receptor alpha, Interleukin 2 receptor alpha chain, IL 2RA, Ly43, p55, TCGFR
Background	CD25 is the alpha chain of the IL-2 receptor. It is a type I transmembrane protein present on activated T cells, activated B cells.