



c-Myb (phospho Ser12) Polyclonal Antibody

Cat No: HR1AP7407

For research use only

Overview

Product Name	c-Myb (phospho Ser12) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	65kDa
GeneID?Human?	MYB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MYB; Transcriptional activator Myb; Proto-oncogene c-Myb
Background	MYB proto-oncogene, transcription factor(MYB) Homo sapiens This gene encodes a protein with three HTH DNA-binding domains that functions as a transcription regulator. This protein plays an essential role in the regulation of hematopoiesis. This gene may be aberrantly expressed or rearranged or undergo translocation in leukemias and lymphomas, and is considered to be an oncogene. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2016],