

GATA-4 (phospho Ser262) Polyclonal Antibody

Cat No: HR1AP6593

For research use only

Overview

Product Name	GATA-4 (phospho Ser262) Polyclonal Antibody
Source	Rabbit
Applications	WB,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GeneID?Human?	GATA4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	GATA4; Transcription factor GATA-4; GATA-binding factor 4
Background	GATA binding protein 4(GATA4) Homo sapiens This gene encodes a member of the GATA family of zinc-finger transcription factors. Members of this family recognize the GATA motif which is present in the promoters of many genes. This protein is thought to regulate genes involved in embryogenesis and in myocardial differentiation and function, and is necessary for normal testicular development. Mutations in this gene have been associated with cardiac septal defects. Additionally, alterations in gene expression have been associated with several cancer types. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2015],