



ZFP598 (phospho Tyr306) Polyclonal Antibody

Cat No: HR1AP9065

For research use only

Overview

Product Name	ZFP598 (phospho Tyr306) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	99kDa
GenID?Human?	ZNF598
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF598; Zinc finger protein 598
Background	<p>zinc finger protein 598(ZNF598) Homo sapiens Zinc-finger proteins bind nucleic acids and play important roles in various cellular functions, including cell proliferation, differentiation, and apoptosis. This protein and Grb10-interacting GYF protein 2 have been identified as components of the mammalian 4EHP (m4EHP) complex. The complex is thought to function as a translation repressor in embryonic development. [provided by RefSeq, Oct 2012],</p>