

Jun D (phospho Ser255) Polyclonal Antibody

Cat No: HR1AP7097

For research use only

Overview

Product Name	Jun D (phospho Ser255) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IP,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	kDa
GeneID?Human?	JUND
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	JUND; Transcription factor jun-D
Background	JunD proto-oncogene, AP-1 transcription factor subunit(JUND) Homo sapiens The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. This protein has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternative translation initiation site usage results in the production of different isoforms (PMID:12105216). [provided by RefSeq, Nov 2013],