



Acetyl Histone H2B (K12) Polyclonal Antibody

Cat No: HR1AP7750

For research use only

Overview

Product Name	Acetyl Histone H2B (K12) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	H2BFS
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	H2BFS; Histone H2B type F-S; Histone H2B.s; H2B/s
Background	H2BFS (H2B Histone Family Member S (Pseudogene)) is a Pseudogene. Diseases associated with H2BFS include endometrial stromal sarcoma. Among its related pathways are Packaging Of Telomere Ends. GO annotations related to this gene include sequence-specific DNA binding and protein heterodimerization activity. An important paralog of this gene is HIST1H2BH.