

ZZZ3 Polyclonal Antibody

Cat No: HR1AP4856

For research use only

Overview

Product Name	ZZZ3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	102kDa
GeneID?Human?	ZZZ3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZZZ3; ZZ-type zinc finger-containing protein 3
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 HTH myb-type DNA-binding domain.,similarity:Contains 1 ZZ-type zinc finger.,