

ZNF695 Polyclonal Antibody

Cat No: HR1AP7980

For research use only

Overview

Product Name	ZNF695 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	ZNF695
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF695; Zinc finger protein 695; Zinc finger protein SBZF3
Background	function:May be involved in transcriptional regulation.,similarity:Contains 1 KRAB domain.,