

ZNF608 Polyclonal Antibody

Cat No: HR1AP4851

For research use only

Overview

Product Name	ZNF608 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	160kDa
GeneID?Human?	ZNF608
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ZNF608; KIAA1281; Zinc finger protein 608; Renal carcinoma antigen NY-REN-36
Background	miscellaneous:Antigen recognized by autologous antibody in patients with renal-cell carcinoma.,PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 C2H2-type zinc finger.,