




## THAP4 Polyclonal Antibody

Cat No: HR1AP4714

For research use only

### Overview

Product Name	THAP4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	63kDa
GenID?Human?	THAP4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	THAP4; CGI-36; PP238; THAP domain-containing protein 4
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Contains 1 THAP-type zinc finger.,



**bioelsa**