

Ribosomal Protein L10L Polyclonal Antibody

Cat No: HR1AP4467

For research use only

Overview

Product Name	Ribosomal Protein L10L Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	20kDa
GeneID?Human?	RPL10L
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RPL10L; 60S ribosomal protein L10-like
Background	ribosomal protein L10 like(RPL10L) Homo sapiens This gene encodes a protein sharing sequence similarity with ribosomal protein L10. It is not currently known whether the encoded protein is a functional ribosomal protein or whether it has evolved a function that is independent of the ribosome. This gene is intronless. [provided by RefSeq, Jul 2008],