




## RECS1 Polyclonal Antibody

Cat No: HR1AP4447

For research use only

### Overview

Product Name	RECS1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	35kDa
GenelD?Human?	TMBIM1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	TMBIM1; LFG3; RECS1; PP1201; PSEC0158; Protein lifeguard 3; Protein RECS1 homolog; Transmembrane BAX inhibitor motif-containing protein 1
Background	similarity:Belongs to the BI1 family.,



# bioelsa