

EFHD1 Mouse mAb

Cat No: HR1AM2196

For research use only

Overview

Product Name	EFHD1 Mouse mAb
Source	Mouse
Applications	WB, IF
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	WB 1:2,000 IF 1:100-200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	27kDa
GeneID?Human?	80303
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	MST133, MSTP133, PP3051, Swiprosin 2, SWS2
Background	EFHD1 is an EF-hand domain-containing protein that displays increased expression during neuronal differentiation.