



## RBM34 Polyclonal Antibody

Cat No: HR1AP6404

For research use only

### Overview

Product Name	RBM34 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GenelD?Human?	RBM34
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	RBM34; KIAA0117; RNA-binding protein 34; RNA-binding motif protein 34
Background	PTM:Phosphorylated upon DNA damage, probably by ATM or ATR.,similarity:Belongs to the RRM RBM34 family.,similarity:Contains 2 RRM (RNA recognition motif) domains.,