

HEI10 Polyclonal Antibody

Cat No: HR1AP3618

For research use only

Overview

Product Name	HEI10 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	32kDa
GeneID?Human?	CCNB1IP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CCNB1IP1; C14orf18; HEI10; E3 ubiquitin-protein ligase CCNB1IP1; Cyclin-B1-interacting protein 1; Human enhancer of invasion 10
Background	cyclin B1 interacting protein 1(CCNB1IP1) Homo sapiens HEI10 is a member of the E3 ubiquitin ligase family and functions in progression of the cell cycle through G(2)/M.[supplied by OMIM, Apr 2004].