



Helios Polyclonal Antibody

Cat No: HR1AP6366

For research use only

Overview

Product Name	Helios Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	56kDa
GeneID?Human?	IKZF2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	IKZF2; HELIOS; ZNFN1A2; Zinc finger protein Helios; Ikaros family zinc finger protein 2
Background	IKAROS family zinc finger 2(IKZF2) Homo sapiens This gene encodes a member of the Ikaros family of zinc-finger proteins. Three members of this protein family (Ikaros, Aiolos and Helios) are hematopoietic-specific transcription factors involved in the regulation of lymphocyte development. This protein forms homo- or hetero-dimers with other Ikaros family members, and is thought to function predominantly in early hematopoietic development. Multiple transcript variants encoding different isoforms have been found for this gene, but the biological validity of some variants has not been determined. [provided by RefSeq, Jul 2008],