



## FoxN4 Polyclonal Antibody

Cat No: HR1AP5809

For research use only

### Overview

Product Name	FoxN4 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GenelD?Human?	FOXN4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXN4; Forkhead box protein N4
Background	forkhead box N4(FOXN4) Homo sapiens Members of the winged-helix/forkhead family of transcription factors, such as FOXN4, are characterized by a 110-amino acid DNA-binding domain that can fold into a variant of the helix-turn-helix motif consisting of 3 alpha helices flanked by 2 large loops or wings. These transcription factors are involved in a variety of biologic processes as key regulators in development and metabolism (Li et al., 2004 [PubMed 15363391]).[supplied by OMIM, Mar 2008],