



FoxD3 Polyclonal Antibody

Cat No: HR1AP3473

For research use only

Overview

Product Name	FoxD3 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	48kDa
GenelD?Human?	FOXD3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXD3; HFH2; Forkhead box protein D3; HNF3/FH transcription factor genesis
Background	forkhead box D3(FOXD3) Homo sapiens This gene belongs to the forkhead family of transcription factors which is characterized by a distinct forkhead domain. Mutations in this gene cause autoimmune susceptibility 1. [provided by RefSeq, Nov 2008],