



CRP1 Polyclonal Antibody

Cat No: HR1AP3165

For research use only

Overview

Product Name	CRP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GenID?Human?	CSRP1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CSRP1; CSRP; CYRP; Cysteine and glycine-rich protein 1; Cysteine-rich protein 1; CRP; CRP1
Background	cysteine and glycine rich protein 1(CSRP1) <i>Homo sapiens</i> This gene encodes a member of the cysteine-rich protein (CRP) family. This gene family includes a group of LIM domain proteins, which may be involved in regulatory processes important for development and cellular differentiation. The LIM/double zinc-finger motif found in this gene product occurs in proteins with critical functions in gene regulation, cell growth, and somatic differentiation. Alternatively spliced transcript variants have been described. [provided by RefSeq, Aug 2010].