

Cytokeratin 6 Mouse mAb Cat No: HR1AM2191

For research use only

Overview

Product Name	Cytokeratin 6 Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:1,000 IHC 1:200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	56-60kDakDa
GeneID?Human?	286887
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	Cytokeratins are proteins of keratin-containing intermediate filaments found in the intracytoplasmic cytoskeleton of epithelial tissue.