

# Jun B (phospho Ser259) Polyclonal Antibody

Cat No: HR1AP2462

For research use only

## Overview

|                       |                                                                                                                                                                                                                                                                                                                                                                                                                             |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | Jun B (phospho Ser259) Polyclonal Antibody                                                                                                                                                                                                                                                                                                                                                                                  |
| Source                | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Applications          | WB,IHC-p,IP,ELISA                                                                                                                                                                                                                                                                                                                                                                                                           |
| Species Reactivity    | Human,Mouse,Rat                                                                                                                                                                                                                                                                                                                                                                                                             |
| Recommended Dilutions |                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Immunogen             |                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Species               | Rabbit                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Storage               | -20°C/1 year                                                                                                                                                                                                                                                                                                                                                                                                                |
| Isotype               |                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Clonality             |                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Concentration         | 1 mg/ml                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Observed band         | 33kDa                                                                                                                                                                                                                                                                                                                                                                                                                       |
| GeneID?Human?         | JUNB                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Human Swiss-Prot No.  |                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Cellular localization |                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Alternative Names     | JUNB; Transcription factor jun-B                                                                                                                                                                                                                                                                                                                                                                                            |
| Background            | function:Transcription factor involved in regulating gene activity following the primary growth factor response. Binds to the DNA sequence 5'-TGA[CG]TCA-3'.induction:By growth factors.,similarity:Belongs to the bZIP family.,similarity:Belongs to the bZIP family. Jun subfamily.,similarity:Contains 1 bZIP domain.,subunit:Binds DNA as an homodimer or as an heterodimer with another member of the jun/fos family., |