



FoxO4 (Acetyl Lys189) Polyclonal Antibody

Cat No: HR1AP9733

For research use only

Overview

Product Name	FoxO4 (Acetyl Lys189) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	FOXO4 AFX AFX1 MLLT7
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	Forkhead box protein O4 (Fork head domain transcription factor AFX1)
Background	forkhead box O4(FOXO4) Homo sapiens This gene encodes a member of the O class of winged helix/forkhead transcription factor family. Proteins encoded by this class are regulated by factors involved in growth and differentiation indicating they play a role in these processes. A translocation involving this gene on chromosome X and the homolog of the Drosophila trithorax gene, encoding a DNA binding protein, located on chromosome 11 is associated with leukemia. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2010],