



ABCB5 Mouse mAb

Cat No: HR1AM2185

For research use only

Overview

Product Name	ABCB5 Mouse mAb
Source	Mouse
Applications	WB, IHC
Species Reactivity	Human
Recommended Dilutions	WB 1:2,000 IHC 1:200-500
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	IgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	90kDa
GeneID?Human?	340273
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	N/A
Background	ATP-binding cassette sub-family B member 5 also known as P-glycoprotein ABCB5 is a plasma membrane-spanning protein that in humans is encoded by the ABCB5 gene. ABCB5 is an ABC transporter and P-glycoprotein family member principally expressed in physiological skin and human malignant melanoma.