

LYAR Polyclonal Antibody

Cat No: HR1AP3844

For research use only

Overview

Product Name	LYAR Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	LYAR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	LYAR; PNAS-5; Cell growth-regulating nucleolar protein
Background	similarity:Contains 2 C2HC-type zinc fingers.,