

Ets-1 (phospho Thr38) Polyclonal Antibody

Cat No: HR1AP2631

For research use only

Overview

Product Name	Ets-1 (phospho Thr38) Polyclonal Antibody
Source	Rabbit
Applications	WB,IHC-p,IF,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	54kDa
GeneID?Human?	ETS1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	ETS1; EWSR2; Protein C-ets-1; p54
Background	ETS proto-oncogene 1, transcription factor(ETS1) Homo sapiens This gene encodes a member of the ETS family of transcription factors, which are defined by the presence of a conserved ETS DNA-binding domain that recognizes the core consensus DNA sequence GGAA/T in target genes. These proteins function either as transcriptional activators or repressors of numerous genes, and are involved in stem cell development, cell senescence and death, and tumorigenesis. Alternatively spliced transcript variants encoding different isoforms have been described for this gene.[provided by RefSeq, Jul 2011],