

HDAC2 (phospho Ser394) Polyclonal Antibody

Cat No: HR1AP2528

For research use only

Overview

Product Name	HDAC2 (phospho Ser394) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat,Monkey
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	55kDa
GeneID?Human?	HDAC2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	HDAC2; Histone deacetylase 2; HD2
Background	histone deacetylase 2(HDAC2) Homo sapiens This gene product belongs to the histone deacetylase family. Histone deacetylases act via the formation of large multiprotein complexes, and are responsible for the deacetylation of lysine residues at the N-terminal regions of core histones (H2A, H2B, H3 and H4). This protein forms transcriptional repressor complexes by associating with many different proteins, including YY1, a mammalian zinc-finger transcription factor. Thus, it plays an important role in transcriptional regulation, cell cycle progression and developmental events. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2010],