



## HNF-3?/?/? (Acetyl Lys264/253/211) Polyclonal Antibody

Cat No: HR1AP2212

For research use only

### Overview

Product Name	HNF-3?/?/? (Acetyl Lys264/253/211) Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GenID?Human?	FOXA1/FOXA2/FOXA3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	FOXA1; HNF3A; TCF3A; Hepatocyte nuclear factor 3-alpha; HNF-3-alpha; HNF-3A; Forkhead box protein A1; Transcription factor 3A; TCF-3A; FOXA2; HNF3B; TCF3B; Hepatocyte nuclear factor 3-beta; HNF-3-beta
Background	forkhead box A1(FOXA1) <i>Homo sapiens</i> This gene encodes a member of the forkhead class of DNA-binding proteins. These hepatocyte nuclear factors are transcriptional activators for liver-specific transcripts such as albumin and transthyretin, and they also interact with chromatin. Similar family members in mice have roles in the regulation of metabolism and in the differentiation of the pancreas and liver. [provided by RefSeq, Jul 2008],