

CD1 Mouse mAb Cat No: HR1AM2179

For research use only

Overview

Product Name	CD1 Mouse mAb
Source	Mouse
Applications	IHC
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	IHC 1:200
Immunogen	
Species	Mouse
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.4. Store at -20°C. Avoid repeated freeze-thaw cycles.
Isotype	lgG1
Clonality	Monoclonal
Concentration	1mg/ml
Observed band	N/AkDa
GeneID?Human?	909
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	CD1A, CD1a molecule, FCB6, HTA1, hTa1 thymocyte antigen, R4, T6
Background	CD1 (cluster of differentiation 1) is a family of glycoproteins expressed on the surface of various human antigen-presenting cells. They are related to the class I MHC molecules, and are involved in the presentation of lipid antigens to T cells.